

Good Partner for the Refractometer !

The 60-C3 circulating constant-temperature bath is an external circulation type constant-temperature water tank provided with a thermo-module. The temperature can be set in the range of 10°C to 60°C and the unit is perfectly suited for connection with a refractometer.

The unit with its compact design can be easily installed anywhere and is easy to operate. The unit is also provided with safety devices such as a low-water cut off device and a thermostat.

Circulating constant temperature bath

60-C3

Cat.No.1921



Specifications

Tank capacity	: 0.5 ℓ
Temperature setting range	: 10 to 60°C(water)
Constant-temperature accuracy	: ±0.2°C or less
Temperature control	: PID control method using a thermo-module
Flow rate	: 4ℓ /min.
Temperature display	: LED (red)
Setting display	: LED (green)
Nozzle	: External diameter of 9 mm, 2 pieces
Temperature sensor	: Platinum thermometer bulb
Safety device	: Overheating preventing device and low-water cut off device
Power supply	: AC 100-240 V, 50/60 Hz
Power consumption	: 430 VA
Dimensions and weight	: 38(W) × 24(D) × 27(H) cm, 11kg
Accessory	: Connection hose (2 m)

All ATAGO products are designed and manufactured in Japan.

HACCP G.M.P. G.L.P.

ATAGO products are acknowledged by HACCP, G.M.P. and G.L.P.

ATAGO CO.,LTD.

Headquarters: 32-10, Honcho, Itabashi-ku, Tokyo 173-0001 Japan
TEL: 81-3-3964-6156 FAX: 81-3-3964-6137
overseas@atago.net <http://www.atago.net/>

ATAGO U.S.A., Inc.

12011 NE First Street, Bldg. C, Suite 110 Bellevue, WA 98005 U.S.A.
TEL: 1-425-637-2107 FAX: 1-425-637-2110
customerservice@atago-usa.com <http://www.atago.net/>

ATAGO INDIA Instruments Pvt. Ltd.

A-501, Mangal Aarambh Business Centre, Near Kora Kendra, Off. S.V. Road, Borivali (West), Mumbai-400 092 India
TEL: 91-22-2833-8038 / 8076 FAX: 91-22-2899-8164
customerservice@atago-india.com



※Specifications and appearance are subject to change without notice.